

Installing The CareEnhance® Review Manager Print Utility

You can install the The CareEnhance® Review Manager Print Utility in two ways:

- ▶ Installing from a portable drive or network shared drive
- ▶ Installing from a CD

This document provides instructions for both.

Note:



These steps must be completed by the local administrator of the workstation.

Installing from a portable drive or network shared drive

1. Unzip the file RMUtilsIE.zip on the drive.
Three files will be displayed: RMUtilsIE.dll, RegisterRMUtilsIE.cmd, and InstallRMUtilsIE.cmd.
2. On the client machine, copy the files RMUtilsIE.dll and RegisterRMUtilsIE.cmd to the local hard drive, for example, C:\McKesson.
3. Execute the RegisterRMUtilsIE.cmd file (if using File Explorer, double-click on the file).
This will unregister the old RMUtilsIE.dll (if one exists) and register the new RMUtilsIE.dll.
4. Test by logging into Review Manager and printing.

Installing from a CD

1. Unzip the file RUtilsIE.zip onto your local drive.
Three files will be displayed: RUtilsIE.dll, RegisterRUtilsIE.cmd, and InstallRUtilsIE.cmd.
2. Burn the files RUtilsIE.dll and InstallRUtilsIE.cmd onto a CD.
3. On the client machine, execute the InstallRUtilsIE.cmd file (if using File Explorer, double-click on the file). This batch file will do the following:
 - ▶ Create the C:\McKesson directory (if one doesn't already exist)
 - ▶ Copy the RUtilsIE.dll file to that directory
 - ▶ Unregister the old RUtilsIE.dll (if one exists)
 - ▶ Register the new RUtilsIE.dll
4. Test by logging into Review Manager and printing.

Note:



If you want to use a directory other than C:\McKesson, you can modify the InstallRUtilsIE.cmd before burning it onto CD. You should be familiar with batch files before making any edits.